

JOBST® Accessories

JOBST Stocking Donner

Product Code: 79668-00

Due to the nature of compression stockings, people may find putting them on to be initially a challenging task. The **JOBST Stocking Donner**, which features a **lightweight frame** and foam grip handles, is designed to aid in the application of all JOBST Medical LegWear.



1 Place your JOBST compression stocking inside the semi-circle frame and pull the top of the stocking down over the semi-circle. The heel of the stocking must face the back of your JOBST® Stocking Donner. Continue to slide the stocking down over the semi-circle until the heel protrudes from the frame, and is centred and even.



2 Whilst using a support for balance, insert your foot into the stocking until your foot is completely on the floor. Helpful Tip: Sit down and hold the heel pocket against the frame as you slide your toes into the stocking.



3 Preferably whilst seated, grasp the padded handles and gently pull the stocking donner upwards whilst moving your leg slightly outwards, until the stocking is above the calf.



4 Push the JOBST Stocking Donner backward and down to free the stocking. Be sure the heel of the stocking is properly placed on the foot. Then adjust the length and smooth out any wrinkles by stroking with the palm of your hands. The knee band in Knee High socks or stockings should stop two finger widths below the crease at the back of the knee.

“IT Stays!”® Roll-on Body Fixative

Product Code: 79274-01

- “IT Stays!”® is a roll-on body fixative used to hold JOBST® LegWear in place
- “IT Stays!”® is pliable, moves with your body and can be washed off easily with water
- “IT Stays!”® does not stain
- Available in a 59ml bottle



JOBST® General Care Instructions

Washing Instructions

- Hand or gentle machine warm wash at 40 degrees celsius (Laundry bag recommended)
- Use mild liquid detergent
- Do not use powder detergent as it clogs elastic fibres
- Rinse well to remove all traces of detergent

For garments with silicone bands

The use of skin care products (such as skin creams, body lotion or powder) as well as various skin types (such as dry or peeling skin) and hair growth, may prevent the silicone band from keeping the thigh high in place. Therefore, moisturisers should be well absorbed before applying stockings.

We recommend the following steps to maximise the grip performance:

1. Turn the garment inside out before washing
2. Before the main wash, thoroughly hand wash the silicone band using a delicate fabric detergent. Concentrate on removing any lotions, hair or foreign material from the silicone.

Drying Instructions

- Dry JOBST® UltraSheer, JOBST® Travel Socks and JOBST® UlcerCARE™ liner either indoors or in the shade, laid flat
- JOBST® forMen Casual, ActiveWear, forMen, Relief and UlcerCARE™ outer stocking can all be tumble dried at low temperature
- Do not wring



Please refer to the product label and / or package insert for full instructions on the safe use of these products.

Always read the label and use only as directed. If symptoms persist, see your doctor / healthcare professional.