Invacare® Softform® Premier Active 2

EN	Mattresses	<i>User Manual</i>
DE	Matratzen	<i>Gebrauchsanweisung</i>
FR	Matelas	<i>Manuel d'utilisation</i>
ES	Colchones	<i>Manual del usuario</i>
IT	Materassi	<i>Manuale d'uso</i>
PT	Colchões	<i>Manual de instruções</i>
SV	Madrasser	<i>Bruksanvisning</i>
NO	Madrasser	<i>Brukerveiledning</i>
DA	Madrasser	<i>Brugsanvisning</i>

This manual must be given to the user of the product.
Before using the product, read this manual and save it for future reference!

Yes, you can.®

- EN:** *Complies with Directive 93/42/EEC concerning medical devices.*
- DE:** *Konform mit Richtlinie 93/42/EWG für Medizinprodukte.*
- FR:** *Conforme à la directive 93/42/ECC relative aux dispositifs médicaux.*
- IT:** *Conforme alla direttiva 93/42/ECC relativa ai prodotti medicali.*
- ES:** *Cumple con la directiva 93/42/EEC relativa a los productos sanitarios.*
- PT:** *Em conformidade com a directiva 93/42/EWG do Conselho relativa aos dispositivos médicos.*
- SV:** *Uppfyller direktiv 93/42/EWG om medicintekniska produkter.*
- NO:** *Overholder retningslinjerne i 93/42/EWG om medisinsk utstyr.*
- DA:** *Overholder direktiv 93/42/EWG om medicinske produkter.*

© Invacare Ltd

Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Invacare Ltd reserves the right to modify existing models at any time. Invacaredealers will be notified of any such modifications. Any use of this publication, or of parts thereof, as well as any reproduction of images, must have the written consent of Invacare Ltd.

Softform® is a registered trademark of Invacare Ltd.

1	General	4	7	After-Use	14
1.1	General information	4	7.1	Storage	14
1.2	Symbols	4	7.2	Re-Use	15
1.3	Symbols Guide	4	7.3	Disposal	15
1.4	Guarantee	4	8	Troubleshooting	15
1.5	Intended use	5	9	Technical Data	16
1.6	Product Description	5	9.1	General Data	16
1.7	Service Life	5	9.2	Pump Unit	16
2	Safety	6			
2.1	Safety information	6			
3	Setup and installation	7			
3.1	Safety information	7			
3.2	Installing the system	8			
4	Usage	8			
4.1	Safety information	8			
4.2	CPR Procedure	9			
4.3	Pump Software Schematic	10			
4.4	Pump Menu Display	11			
5	Transport	13			
5.1	Safety information	13			
5.2	Transport Mode	13			
6	Maintenance	13			
6.1	Inspection	13			
6.2	Cleaning and Care	13			

I General

I.1 General information

Essential nursing care is pivotal in pressure ulcer prevention. The Softform® Premier Active 2 mattress will positively contribute to the outcome of a pressure ulcer prevention care plan.

Education, clinical judgement and action based planning based on vulnerability are fundamental factors in prevention of pressure ulcers.

A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

Please heed all the notes, particularly the safety information, and act accordingly.

Should you require more information on the Invacare® Softform® Premier Active 2 mattress please refer to the Invacare® Softform® Brochure and contact details on the back page of this User Manual.

I.2 Symbols

Warnings

In this User Manual warnings are indicated by symbols. The warning symbols are accompanied by a heading that indicates the severity of the danger.

WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

I.3 Symbols Guide

	Do not pierce or cut		Line dry
	Do not dry clean		Declaration of conformity
	Do not put near flame		Recommended 80 °C
	User weight limit of 247.6 kg		Do not iron
	Do not bleach		Tumble dry low heat

I.4 Guarantee

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business. Guarantee claims can only be made through the relevant specialist dealer.

Standard Invacare® Terms

This is to certify that your Softform® Premier Active 2 Mattress is warranted by Invacare® Ltd for a period stated in the Table "Technical Data" of this user guide.

The Warranty of your Invacare® Softform® product is valid from time of shipping.

If a defect or fault is discovered the Invacare® dealer or Local Business Development Manager from whom the appliance was obtained must be notified immediately.

The manufacturer will not accept responsibility for damage caused by misuse or non-observance of the instructions set out in this user guide.

During the period of the warranty any products that have become defective due to faulty workmanship or materials will be renewed without charge.

The warranty will be forfeited should any unauthorized alteration be made to the equipment.

Both warranty and fire retardancy Certification will become null and void if non-Invacare® spares are used on any Invacare® Softform® Mattress products.

The purchaser's statutory rights under the Consumer Protection Act are not affected.

Quality and Flame Retardancy

Quality is fundamental to the company's operation, working within the disciplines of ISO 9001 and ISO 13485.

The Softform® Premier Active 2 Mattress features the CE mark, in compliance with the Medical Device Directive.

The foam and covers used to manufacture the Invacare® Softform® range of mattresses are independently tested and certified in accordance with EN 597-1 and EN 597-2 and BS 7177 Crib 5.

The pump unit is tested to EC Directive 89/336/EEC, EN 55022 and EN 61001-1, 2. Manufactured to comply with EN 60601-1.

For further information please contact Invacare in your country (addresses see back page of this manual).

1.5 Intended use

This pressure redistribution mattress and pump is intended to be used in conjunction with an appropriately sized bed frame, as part of an overall pressure ulcer prevention program of care.

It can be used safely in static mode (deflated) for static pressure redistribution, or in dynamic mode (inflated) should an alternating pressure support surface be required.

This product has been designed to deliver effective pressure reduction to users, when the product is in normal use which is defined by Invacare® Ltd as when the support surface is covered with a cotton, cotton combination or linen bed sheet, and any one of these would be the only item deployed between the support surface and the user.

1.6 Product Description

The Softform® Premier Active 2 mattress system acts as a static pressure reducing support/mattress for patients at high risk that can, by facilitating the air pump, introduce effective alternating pressure if the patients condition requires alternating pressure therapy.

The water-resistant cover provides a vapour-permeable, multi stretch surface, to promote patient comfort and to maximise the effectiveness of the foam core.

The complete system comprises of:

A mattress insert comprised of alternate air cells that connect to a twin air connector hose.

A Micro-processor controlled air pump unit.

A 5 meter power lead.

1.7 Service Life

We estimate a life expectancy of five years for these products, provided they are used in strict accordance with the intended use as set out in this document and all maintenance and service requirements are met. The estimated life expectancy can be exceeded if the product is carefully used and properly maintained, and provided technical and scientific advances do not result in technical limitations. The life expectancy can also be considerably reduced by extreme or incorrect usage.

The fact that we estimate a life expectancy for these products does not constitute an additional warranty.

2 Safety

2.1 Safety information

WARNING

▶ Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals, service manuals or instruction sheets supplied with this product or optional equipment. Invacare product manuals are available at www.invacare.co.uk or your local dealer. If you are unable to understand the warnings, cautions or instructions, please contact a healthcare professional, dealer or technical personnel before attempting to use this equipment – otherwise, injury or damage may occur.

WARNING

Invacare products are specifically designed and manufactured for use in conjunction with Invacare accessories. Accessories designed by other manufacturers have not been tested by Invacare and are not recommended for use with Invacare products.

The introduction of certain third party products between the mattress surface and the user may reduce or impede the clinical effectiveness of this product.

'Third party products' may include, but are not limited to items including under blankets, plastic sheets and sheepskins, etc. Heated over blankets must only be used in consultation with a suitably qualified healthcare professional, as an increase in temperature can increase the risk of developing pressure ulcers.

WARNING

Bed sheets must be loosely fitted, with creases smoothed out. Care must always be taken to ensure that the support surface in contact with the user is kept free from crumbs and other food debris, and that drip cables, stents, and other foreign objects do not become entrapped between the user and the pressure reducing surface of the mattress, as this may result in the development of pressure ulcers.

WARNING Risk of fire or explosion!

A cigarette can burn a hole in the bed surface and cause damage to the mattress. Also, patient clothing, bed sheets, etc, may be combustible and cause a fire. Failure to observe this warning can result in a severe fire, property damage and cause physical injury or death.

There is an explosion risk if used with flammable anesthetics.

There is a possible fire hazard when used with oxygen administering equipment other than nasal mask or half bed tent type.

- ▶ Do not smoke while using this device.
- ▶ An oxygen tent may not extend below mattress support level.

IMPORTANT

The information contained in this document is subject to change without notice.

- ▶ Check all parts for shipping damage and test before using.
- ▶ In case of damage, do not use.
- ▶ Contact Invacare/Carrier for further instructions.

3 Setup and installation

3.1 Safety information

**DANGER** Electrical shock hazard!

- ▶ Do not remove pump cover.
- ▶ Refer to qualified service personnel.
- ▶ Before performing any maintenance to the power unit, disconnect the power lead from the wall outlet.
- ▶ Do not insert items into any openings of the control unit. Doing so may cause fire or electric shock by shorting the internal components.
- ▶ The control unit must be kept away from all heat sources and radiators during operation.
- ▶ Connect the equipment to properly grounded three prong wall outlet using the 5 meter power lead provided with the product.
- ▶ Grounding reliability depends upon a properly grounded receptacle (3-prong).

**WARNING** Risk of entrapment!

Patient entrapment with the bed side rails may cause injury or death.

Proper patient assessment and monitoring, and proper maintenance and use of equipment is required to reduce the risk of entrapment. Variations in bed rail dimensions, and mattress thickness, size and density could increase the risk of entrapment.

- ▶ Mattress must fit bed frame and side rails snugly to prevent patient entrapment. Follow the bed manufacturer's instructions.
- ▶ After any adjustments, repair or service and before use, make sure all attaching hardware is tightened securely. Rails with dimensions different from the original equipment supplied or specified by the bed manufacturer may not be interchangeable and may result in entrapment or other injury.

**WARNING**

The Softform® Premier Active 2 is recommended to be installed on medical bed frames with bed sides or assist rails. It is preferred that the rails to be in the raised position whenever the patient is on the bed. Healthcare professionals assigned to each case should make the final determination whether side assist rails are warranted after assessing patient risk of entrapment.

Controls on the footboard may be obstructed by the power unit on a few bed frames. It may be necessary to relocate the power unit.

- ▶ Before placing the patient on the bed, check that air hoses and power cord are clear of moving bed components.
- ▶ Operate all bed frame motorized functions through their full range of motion to be certain that there is no pulling, interference or pinching.

3.2 Installing the system

	<ul style="list-style-type: none"> ▶ Place the pump unit at the end of the bed (preferably hung by the built in pump hooks).
	<ul style="list-style-type: none"> ▶ Detach the Softform® Premier Active 2 twin hose from the end of the mattress. ▶ Connect it to the pump. (An automatic click will signify a secure connection).
	<ul style="list-style-type: none"> ▶ Connect the mains power lead to the pump unit and a suitable power supply.
	<ul style="list-style-type: none"> ▶ Switch on mains power. ▶ Switch on pump switch. <p>System will pressurise (see chart 4.3, page 10) shown by a green light and three audible beeps.</p> <ul style="list-style-type: none"> ▶ Correct pressure shown by a green light and two audible beeps.

4 Usage

4.1 Safety information

- ▶ Remove all packaging before use.
- ▶ Place the mattress directly on the frame of the bed.

WARNING

It is very important for the patient to reposition themselves, or to be repositioned, on a regular basis. This must be based on the clinical judgement of a qualified healthcare professional. This relieves pressure which helps prevent both tissue compression and potential ulcer formation.

- ▶ Always consult a qualified healthcare professional before using the Softform® Premier Active 2.
- ▶ Monitor the patient frequently.

CAUTION

- ▶ The printed side of the mattress cover must always be facing upwards.
- ▶ Make sure that the distance between the surface of the mattress and the top of the side rail should be at least 220 mm.

IMPORTANT

Medical equipment including infusion pumps and monitors should be attached to appropriate bed accessories. In domestic settings common causes of damage include cigarette burns and the claws of pets that puncture covers, allowing fluid ingress and staining.

! **IMPORTANT**

- To prevent accidental cover damage, do not place hypodermic needles, venflons, scalpels or other similarly sharp objects onto the mattress. Ensure that all venflons are taped down correctly with no sharp edges exposed.
- When using bridging boards or other patient transfer aids, care should be taken not to damage the mattress cover. All transfer aids should be checked for any sharp edges or burrs before use as these can damage the mattress cover.
- Make sure that the mattresses are not jammed or damaged by sharp edges when used on beds with an adjustable frame.
- When using the Softform® Premier Active 2 on a profiling bed ensure that the knee break is used before the backrest.

4.2 CPR Procedure

The Softform® Premier Active 2 mattress has been fully tested to comply with the current CPR standard of 4 – 5 cm compression depth. This was achieved at all stages of inflation/deflation.


A report on these findings conducted by the Resuscitation Unit at the University Hospital of Wales, Cardiff is available upon request.

- Hold down Red CPR button.
- Pull hose connector firmly away from the pump.
- Switch off the pump.


Mattress will start to deflate.


- When CPR is complete reactivate the system following chapter 3.2, page 8.

4.3 Pump Software Schematic

4.4 Pump Menu Display

		Task	Action	Display
1		Power up.	▶ Switch pump on.	Three audible beeps. LCD illuminates. Green light illuminates. “System setup” displayed.
				After 10 minutes “System OK” displayed.
2		Menu option	▶ Cannot enter menu.	Wait until set up is complete “System OK”.
3		Disconnection of CPR or Disconnection of hose in mattress.	▶ Low Pressure alarm.*	“Low Pressure” displayed with audible beeps.
4		Enter menu mode. (3 second time limit between each function).	▶ Press arrow Up.	“System OK” with “1” underneath.
			▶ Press arrow down.	“System OK” with “2” underneath.
			▶ Press OK.	“Pressure” displayed. Pump display will revert back to “System OK” if no selection made within 15 sec.

		Task	Action	Display
5		Pressure adjustment.	<ul style="list-style-type: none"> ▶ Press OK button to select Pressure. 	Default set at “60mmHg”.
			<ul style="list-style-type: none"> ▶ Arrow up to selected pressure. <p>or</p> <ul style="list-style-type: none"> ▶ Arrow down to selected pressure. 	60 ↑ 70 ↑ 80mmHg. 60 ↓ 50 ↓ 40mmHg.
			<ul style="list-style-type: none"> ▶ Press OK button. 	Pressure displayed under “System OK”.
6		Follow sequence 4 to enter Language mode.	<ul style="list-style-type: none"> ▶ Toggle from Pressure using ↑ or ↓ button. ▶ Press OK button. 	“Language” displayed.
7		Language selection using sequence 6.	<ul style="list-style-type: none"> ▶ Toggle from Language using ↑ or ↓ button to desired language. ▶ Press OK button. 	“English” default. Available languages, see chapter 4.3, page 10.

* Low Pressure Alarm will only function after the system has gone through full cycle set up (System OK).

5 Transport

5.1 Safety information

! IMPORTANT

- Take care when handling mattresses to ensure no damage to the cover. It is recommended that two people lift/carry mattresses.
- Avoid contact with jewellery, nails, abrasive surfaces etc.
- Do not drag mattresses.
- Avoid contact with wall, door frames, door catches or locks etc.
- Do not transport in roll cages unless completely protected from the sharp edges of the cage.

5.2 Transport Mode

If it is necessary to move the bed or mattress simply.

- Turn off power supply.
- Disconnect pump lead (if necessary the air hose).
- When system is ready to reactivate following chapter 3.2, page 8.

Air supply hose should be stored by attaching fastener at foot of the mattress.

6 Maintenance

6.1 Inspection

Check mattresses (foam, air insert and cover) for strike-through (this may include fluid ingress, stains, rips or damage) after the release of each patient, after ending of the period of use or on a minimum monthly basis (depending on which occurs first) by a suitably qualified and competent person.

Check mattresses

- Unzip the cover completely.
- Check for any staining on the white underside of the cover.
- Check for any staining on the interior foam.

- Replace any stained items and dispose of as per local authority procedure.

6.2 Cleaning and Care

! IMPORTANT

All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean. For further information on decontamination in Healthcare Environments, please refer to 'The National Institute for 'Clinical Excellence' guidelines on Infection Control www.nice.org.uk/CG002 and your local infection control policy.

Cleaning covers

(Removal of contaminants such as dust and organic matter)

- Remove all covers for laundering.
- Launder the covers with the recommended temperature between 60 and 80 °C using a diluted detergent solution (Instructions on label).

! IMPORTANT

Washing at higher temperatures may cause shrinkage.

Drying covers

- Hang mattress covers from a line or bar and drip dry in a clean indoor environment.

or

- Tumble dry on a low heat setting.

! IMPORTANT

- Tumble dry setting must not exceed 40 °C.
- Do not tumble dry for longer than 10 minutes.
- Dry covers thoroughly before re-fitting to the mattress.

Disinfecting covers

(Reducing the number of microorganisms)

Please contact your hygiene specialist in the event of contamination.

! IMPORTANT

Ensure that any residual detergent has been removed prior to disinfection.

Light soilage

- Wipe down the cover with a 0.1% Chlorine Solution (1,000 ppm).
- Rinse the cover with clean water using a single use nonabrasive cloth.
- Dry the cover thoroughly.

Heavy soilage

Where the mattress is badly soiled, we recommend cleaning with a dilute cleaning solution at 80 °C in the washing machine.

! IMPORTANT

1% Chlorine Solution used on a regular basis can diminish the life of the cover if not rinsed and dried properly.

- Do not use granules.

- Clean up all spillages of bodily fluids i.e. blood, urine, faeces, sputum, wound exudate and all other bodily secretions as soon as possible using a 1% Chlorine Solution (10,000 ppm).

Large spillages of blood should be absorbed and removed with paper towels first and followed by as above.

- Rinse the cover with clean water using a single use nonabrasive cloth, and thoroughly dried.

! IMPORTANT

Polyurethane coated fabrics can absorb liquids for short periods causing a temporary change to the polyurethane characteristics. The mattress cover swells temporarily and is more vulnerable to physical damage for a period after it is completely surface dried, by which time it will revert to its previous state.

- Dry the cover thoroughly.

WARNING

- Remove contaminated foams from use.

CAUTION

- Keep clear of open heat sources.

IMPORTANT

- Do not use phenols, alcohols, bleaches, or other abrasive materials.

Change covers

- Unzip the cover and remove it carefully from the foam core.

IMPORTANT

- Ensure that the corners of the foam core are positioned correctly into the corners of the cover. Ensure that the profiled side of the foam is facing uppermost when packed into its cover.

- Then close the zipper.

7 After-Use

7.1 Storage

IMPORTANT

- Store mattresses in a dry environment.
- Store mattresses within a protective cover.
- Store items in a flat state on clean, dry, off-flooring free from sharp edges to avoid any possible damage.
- Never store other items on top of a mattress.
- Do not store mattresses next to radiators or other heating devices.

7.2 Re-Use

A cleaning record must be kept as part of the cleaning system.

The product is suitable for repeated use. The number of times it can be used depends on how often and in which way the product is used.

▸ Before reuse, clean the product thoroughly, see chapter 6.2, page 13.

7.3 Disposal

The disposal and recycling of used devices and packaging must comply with the applicable legal regulation.

8 Troubleshooting

Problem / Alarm	Cause	Solution
Mattress not inflating (not alternating properly).	Mattress hose disconnected.	Connect hose connectors, lock them in place.
	Air hose kinked or split.	Unkink or replace split hose.
	Major leak in air cell.	Replace leaking air cell.
	Has power and fuse is good, control unit does not come on.	Send control unit back to factory for repair.
	Not alternating, rotor malfunction.	Send control unit back to factory for repair.
	No air (pump malfunction).	Send control unit back to factory for repair.
No power.	Control unit off.	Check power source, turn unit on.
	Power cord disconnected.	Connect power cord and the power source is on.
	No power in the power source.	Check the power source has power and turn it on.
	Power outage.	Wait till the power source has power.
	Control unit switched on but no power to unit, fuse blown.	Consult engineer, refer to service manual.

9 Technical Data

9.1 General Data

Product	Warranty ¹⁾	Firetesting	Grade ref & Colour	Nominal density range [kg/m ³]	Nominal hardness range [N]	Maximum user weight [kg]	Weight of product [kg] ²⁾
SOFTFORM® PREMIER ACTIVE 2	Cover: 4 years Foam: 8 years	EN 597-1, 2 BS 7177: Crib 5	RX 39/120 Pink RX 39/200 Blue	38 – 40 38 – 40	105 – 135 180 – 200	247.6	14

¹⁾ Warranty is against manufacturing defects.

²⁾ Based on the weight of a standard size mattress. This can change if different sizes are ordered.

9.2 Pump Unit

Main Supply	220 / 240 V AC, 50 Hz
Rated Input Current	1 A
Supply Fuse	1 A
Noise Level	< 32 dB
Classification	Class 2 Type BF
Cycle Time	10 min
Size	237 mm x 205 mm x 80 mm
Weight	1.75 kg
Air Flow	4 l/min
Operating Pressure	60 mmHg (8 kPa)
Power	10 W
Warranty on pump	24 months
Warranty on air cell insert	24 months

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Invacare AB, Fagerstagatan 9,
Invacare UK Operations Limited