


Max Mobility

..... live life to the max

DELTA T14 FOLDING COMMODE



Features:

- steel frame & moulded armrests
- slip resistant rubber feet / tips
- height adjustable steel legs
- includes pan, lid & splash-guard
- no tools required for assembly
- user weight capacity of 100kg

DELTA T14 FOLDING COMMODE	Metric	Imperial
Approximate Total Product / Unit Weight	6.4 kg	14.1 lbs
Overall Width (product or frame when opened)	570 mm	22.4 inch
Overall Depth (product or frame when opened)	540 mm	21.3 inch
Width / Depth (product or frame when closed)	120 mm	4.7 inch
Overall Height (or height range if adjustable)	615 - 765 mm	24.2 - 30.1 inch
Armrest Width (outside dimension)	560 mm	22.0 inch
Armrest Width (inside dimension)	455 mm	17.9 inch
Armrest Height (from top of seat)	200 mm	7.9 inch
Back / Back-Rest Height (from top of seat)	200 mm	7.9 inch
Seat Height (from top of seat to floor)	415 - 565 mm	16.3 - 22.2 inch
Seat Width (actual seat, not available space)	350 mm	13.8 inch
Seat Depth (actual seat, not available space)	420 mm	16.5 inch
Wheels (dimensions of front)	N / A	N / A
Wheels (dimensions of rear)	N / A	N / A
Recommended Maximum User Weight	100 kg	220.5 lbs

Cleaning & Maintenance: This product may be cleaned with a soft cloth. The use of water or cleaning products is not recommended as they may damage the product or cause rust. Any wheels, brakes, adjustment pins/knobs or any other moving parts should be checked regularly to ensure they are working correctly.

Warranty & Consumer Guarantee: This product has a 12 month warranty from date of purchase. Valid proof of purchase is required for all warranty claims. The warranty covers manufacturing faults but does not cover change of mind or damage, fault or failure due to external causes, including accident, abuse, misuse, or normal wear and tear. Any unauthorised repairs, alterations or modifications, or any attempt to open or dismantle the product will render the warranty void. The defective product, together with valid proof of purchase, should be returned to the place of purchase or shipped directly to the supplier. If shipped directly to the supplier, you must prepay the shipping and handling charges. The supplier will determine at no cost to you the cause of the fault (if any) within a reasonable time frame. Any repaired or replaced product, original product determined to be not faulty, or refund (if a major failure is determined) will be sent to the place of purchase you returned the product to for collection, or shipped to you freight prepaid if you shipped the product directly to the supplier. This product comes with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the product repaired or replaced if the product fails to be of acceptable quality and the failure does not amount to a major failure. A product presented for repair may be replaced by a refurbished product of the same type rather than being repaired.

Imported by: Home Medical Products Australia Pty Ltd
For: Max Mobility Pty Ltd
461 Whitehorse Road, Balwyn VIC 3103
Phone: (03) 9830 0305

