


ACTIMOVE®

Lower Leg Walker

APPLICATION



ORDERING INFORMATION

Ref-No.	Size	UK	US	EU	
Actimove®					
Standard Walker	Air Walker	Shoe Size			
76270-35	76272-35	XS	<3.5 <3	<4.5 <6	<36.5 <36
76270-36	76272-36	S	4-6.5 3.5-5.5	4.5-7 6-8	36-39.5 36.5-38.5
76270-37	76272-37	M	7-9.5 6-8.5	7.5-10 8.5-11	40-44 39-43
76270-38	76272-38	L	10-12 9-11	10.5-12.5 11.5-13.5	44.5-46.5 44-46
76270-39	76272-39	XL	>12 >11	>12.5 >13.5	>46.5 >46



Actimove®

Feel good. Move better.



EFFECTIVE
IMMOBILISATION,
ENHANCED
PATIENT COMFORT
THE ACTIMOVE® WALKER RANGE

BSN medical Australia & New Zealand

L6, 211 Wellington Road,
Mulgrave, Victoria, Australia 3170

www.bsnmedical.com.au
www.bsnmedical.co.nz

BSN0451

05/2018

Customer Service Australia
T 1300 276 633 F 1300 998 830

Customer Service New Zealand
T 0508 276 111 F 0508 998 830

BSNmedical

an Essity company

THERAPIES. HAND IN HAND.

BSNmedical

an Essity company

ACTIMOVE® LOWER LEG WALKER – ENHANCED FEATURES FOR MAXIMUM CLINICAL BENEFIT

ACTIMOVE® STANDARD AND AIR WALKER

Fractures and soft tissue injuries of the foot and ankle can be treated with a variety of conservative methods. For indications requiring total immobilisation, the **Delta-Cast®** Solutions range offers a comprehensive selection of casting and immobilisation materials. Where a removable device is indicated, the **Actimove® Walker** will stabilise the injury whilst allowing the freedom of normal ambulation and removability.

ACTIMOVE® WALKER

The **Actimove® Walker** range offers impressive features. The high-quality rigid outer skeleton ensures high protection of the foot, ankle and lower leg. Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient. The rocker bottom of the **Actimove® Walker** facilitates a natural physiological gait when walking.



ACTIMOVE® AIR WALKER

The **Actimove® Air Walker** has a liner which can be inflated for oedema control. In just a few pumps, the inflated liner can reduce the void around the foot and ankle. This improves intimacy of fit, reduces unwanted movement, and allows for optimal patient comfort and wearing tolerance.



ADVANCED FEATURES FOR SAFETY AND PATIENT CARE



..... Circumferential strapping
For increased stability.



..... Breathable liner
For optimal skin climate and high wearing comfort.



..... Rocker bottom and shock absorbing sole
Securely simulates a natural gait.

INTEGRATED PUMP SYSTEM OF ACTIMOVE® AIR WALKER

Integrated pump and valve system
Designed to easily adjust desired compression level.



Indications

- Fractures on the foot and ankle
- Acute ankle sprains of the lower leg
- Stress fractures of the lower leg
- Post-operative and rehabilitation use
- Lower extremity soft tissue injuries
- Severe ankle and foot sprains