

Lower Extremity

JOBST® FarrowWrap® Foot - Suitable For 24-Hour Wear



LITE Footpiece

LITE footpiece provides compression for mild to moderate oedema.

- 20-30 mmHg compression

STRONG Footpiece

STRONG Footpiece provides compression for moderate to severe oedema. The long ankle bands can be wrapped around the ankle or over the forefoot to provide more compression where needed.

- 30-40 mmHg compression

CLASSIC Footpiece

CLASSIC footpiece is clinically proven in managing oedema. Durable construction.

- 30-40 mmHg compression

Customising garments

On triangular shaped legs we recommend **JOBST® Farrow GarmentGrip™** be worn over our legpiece to help hold the garment in place.

Lower Extremity

JOBST® Farrow Toe Cap


Ready to wear glove which can be worn on either right or left foot.

Fine seams, with ultra thin fabric, allows air to pass to enhance patient comfort.

The toes can be easily trimmed with scissors to exactly the required length without unravelling.

Measurement Table

Size	cA
XS	20-22 cm
S	22-24 cm
M	24-26 cm
L	26-28 cm



Ordering Information

Product Code	Size	Compression	Colour
72947-00	XS	20 - 30 mmHg	Beige
72947-01	S	20 - 30 mmHg	Beige
72947-02	M	20 - 30 mmHg	Beige
72947-03	L	20 - 30 mmHg	Beige

