

Lower Extremity

JOBST® FarrowWrap® Kneepiece



Adjustable kneepiece, easy to apply & comfortable to wear.

Can be worn in combination with JOBST® FarrowWrap® Legpieces and Thighpieces, or as a single unit (e.g. in combination with an Elvarex knee high stocking)

Available in two fabrics: LITE and STRONG



How to measure

JOBST® FarrowWrap® Kneepiece sizing chart

	Xsmall	Small	Medium	Large	X Large
5cm above patella	47 - 52 cm	53 - 59 cm	60 - 66 cm	67 - 73 cm	74 - 80 cm
Calf	36 - 43 cm	42 - 50 cm	48 - 58 cm	53 - 63 cm	58 - 68 cm

Ordering Information

	Ref - No.	Size	Colour
JOBST® FarrowWrap® Lite (20-30mmHg) Kneepiece			
	73367-00	XSMALL	TAN
	73367-01	SMALL	TAN
	73367-02	MEDIUM	TAN
	73367-03	LARGE	TAN
	73367-04	XLARGE	TAN
JOBST® FarrowWrap® Strong (30-40mmHg) Kneepiece			
	73368-00	XSMALL	TAN
	73368-01	SMALL	TAN
	73368-02	MEDIUM	TAN
	73368-03	LARGE	TAN
	73368-04	XLARGE	TAN

JOBST FarrowWrap® RTW Footpiece size chart

		XSmall	Small	Medium	Large	XLarge
Ready-to-wear pieces						
Footpiece						
A'	Mid-foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm
X	Regular	15-17 cm	16-19 cm	18-20 cm	19-22 cm	21-24 cm
X	Long	18-20 cm	20-22 cm	21-23 cm	23-25 cm	25-27 cm

The length is measured from the most posterior aspect of the heel to the base of the 5th toe.

JOBST FarrowWrap® RTW Legpiece size chart

		XSmall	Small	Medium	Large	XLarge
Ready-to-wear pieces						
Legpiece						
C	Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
B	Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
A-D	Regular	35-39 cm	37-41 cm	39-43 cm	41-45 cm	42-46 cm
A-D	Tall	40-44 cm	42-46 cm	44-48 cm	46-50 cm	47-51 cm

Length A-D = Floor to 2 finger widths below crease of knee. Measure posteriorly FOLLOWING contour of leg.