

JOBST® Bella™ Lite
Ready-to-Wear Upper Limb Garment




**Ready-to-Wear
Combination
Sleeve**



The management of oedema with compression therapy can be challenging as limbs do not follow a normal anatomical shape.

The effective management of oedema is highly dependent on a well-fitting garment that contours to the individual patient’s anatomical differences/shape. Custom garments provide excellent fit and are therapeutically effective. **JOBST® Elvarex® Custom** garments is recognised as a preferred option for these patients.

There are however situations where the patient’s symptoms severity and/or progression will allow flexibility in treatment options. The **JOBST®** range provides a quality garment option for all of these patients.



Early intervention/prevention	Mild lymphoedema	Moderate lymphoedema	Severe lymphoedema
Bella™ Lite	Bella™ Lite – Good	Bella™ Strong – Good	Elvarex® Custom Made
	Bella™ Strong – Better	Elvarex® RTW – Better	
		Elvarex® Custom Made – Best	



Introducing Bella™ Lite
Combination Sleeve,
providing the comfort of
an **all in one** garment
for your patients.

JOBST® Bella™ Lite Combined Sleeve

Bella™ Lite Combination Sleeve provides
the comfort and effectiveness of the
Bella™ Lite Armsleeve and Gauntlet in an
all in one garment



**Measurement
made easy with
3 simple steps**



**Affordable
treatment option
for your patients**



**Ideal as an early
intervention
treatment option**



Early intervention/prevention

Bella™ Lite™

JOBST® Bella™ Lite

Ready-to-Wear Upper Limb Garment

JOBST® Ready-to-Wear Bella™ Lite helps maintain the limb size of lymphoedema patients. These comfortable garments provide graduated compression to assist in patient management following trauma, surgery or radiation therapy.

JOBST® Bella™ Lite

- Soft and air permeable
- Lightweight, knit construction
- Fits a wide range of arm and hand shapes
- Easy to don, comfortable to wear
- Armsleeve and Gauntlet 20-30mmHg

2 lengths for optimal fit



Even more improvement:

All sizes are available in 2 lengths:

Regular and Long



The JOBST® Bella™ Lite range

Indications:

- Oedema
- Lymphoedema

The armsleeve can also be used for:




- Phlebitis
- Post-thrombotic syndrome
- Vascular malformations


Contraindications:

- Arterial insufficiency
- Cutaneous infection
- Dermatitis in the acute phase
- Wet dermatosis
- Conditions where increased venous and/or lymphatic return is not desired

Our new lengths are based on measuring data from over 36,000 patients.

JOBST® Bella™ Lite combination sleeve measurement information/sizing charts

Size	Code	 Palm	 Mid-Lower	 Mid-Upper
S	77680-00 (reg) / 77680-03 (long)	15cm - 19cm	18cm - 23cm	22cm - 40cm
M	77680-01 (reg) / 77680-04 (long)	19cm - 22cm	23cm - 29cm	27cm - 45cm
L	77680-02 (reg) / 77680-05 (long)	22cm - 26cm	29cm - 34 cm	33cm - 51cm

Lengths	 C	 G
Regular (up to 43cm)		Long (over 43cm)

DISTRIBUTED BY:

Smith & Nephew Pty Ltd Healthcare Division

Australia

T 61 3 8540 6777 F 61 3 9544 5086
www.smith-nephew.com/australia

Customer Service

T 13 13 60 F 1800 671 000

New Zealand

T 64 9 828 4059 F 64 9 820 2867
www.smith-nephew.com/new-zealand

Customer Service

T 0800 807 663 F 0800 263 222



©Registered Trademark
SN11864 (02/15)