

The Zoid Cushion



1. Air-Foam Floatation

Functional areas of the cushion are created by three types of foam; soft foam for the decubitus-sensitive area of the ITs; medium foam for main seat support area; and firm foam for the perimeter.

2. Positioning Supports

Extra-firm, bevelled foam increases lateral stability and sitting balance, while a large ischial pan protects the pelvis in both anterior and posterior pelvic tilt positions. There is no medial thigh separator, therefore lower extremities can come into midline for better fit in a tapered chair.

3. Valve for Adjusting Air Volume

Zoid PSV is supplied with Varilite's revolutionary Pressure Setting Valve.

4. Cover

Stylish and functional. Available in mesh and incontinence. Zoid PSV covers meet ISO 7176-16 ignition resistance standards for upholstered wheelchair components.



Zoid

	12x8" (30x35cm)	12x12" (35x35cm)	12x14" (35x40cm)	14x14" (35x45cm)	14x16" (35x50cm)	14x18" (38x38cm)	15x15" (38x43cm)	15x17" (40x35cm)	16x16" (40x40cm)
w/ Mesh Cover	32800	32200	32400	34400	34600	34800	35500	35700	36600
w/ Incontinence Cover	32810	32210	32410	34410	34610	34810	35510	35710	36610
	16x18" (40x45cm)	16x20" (40x50cm)	17x17" (43x43cm)	18x16" (45x40cm)	18x18" (45x45cm)	18x20" (45x50cm)	20x16" (50x40cm)	20x18" (50x45cm)	20x20" (50x50cm)
w/ Mesh Cover	36800	36000	37700	38300	38800	38000	30600	30800	30000
w/ Incontinence Cover	36810	36010	37710	38310	38810	38010	30610	30810	30010

Zoid Spare Covers

	12x8" (30x35cm)	12x12" (35x35cm)	12x14" (35x40cm)	14x14" (35x45cm)	14x16" (35x50cm)	14x18" (38x38cm)	15x15" (38x43cm)	15x17" (40x35cm)	16x16" (40x40cm)
Mesh Cover	04651	04652	04653	04654	04655	04697	04698	04699	